

Grip Ribbon

FOR SUPPORT RAILS AND HANDLES

S

length 45 cm, width max. 34 cm

L

length 45 cm, width max. 41 cm



Increased friction on support rails and handles

- When grip strength weakens, a slight increase in friction enhances user autonomy and independence.
- Can be machine washed at 60 degrees Celsius with mild detergent, rinsing is recommended.
 - Use at room temperature.
- Made of synthetic rubber, does not contain latex or silicone.
 - Dispose as general waste.

Helps to increase user's functional capacity and independence