

Grip Ribbon

STEADYHAND

S

length 45 cm, width max. 34 cm

L

length 45 cm, width max. 41 cm



Assistive device for maintaining the position of a neurological patient

- Increased friction enhances grip – useful e.g. while eating and moving around with mobility aids.
- Breathable mesh tube that adapts to the shape of the hand.
 - The mesh structure of the product does not unravel.
- Can be machine washed at 60 degrees Celsius with mild detergent, rinsing is recommended.
 - Use at room temperature.
 - Dispose as general waste.

Cut a hole for the thumb in the tubular ribbon to suit its intended use